



BNX Converting LLC. Solutions Brief

Supply of N95 Respirator Face Masks
Area of Interest (Aol)

2022

CAGE Code: 8T7Z3

DUNS number: 117714784

NAICS: 339113, 423450, 561990, 624230

Socio-Economic Status: Minority Owned Business

CONTACT

Adam Ravat

adam.ravat@bnx.com

315-790-0466

BNX Converting LLC.

16727 Park Row,
Houston, Texas,
77084

Company Overview

BNX Converting was launched in 2020 by Amcrest, which has been involved in the manufacture and distribution of security cameras since 2010.

BNX Converting is a limited liability company that was founded in Houston, Texas in order to support the public response to the COVID-19 outbreak.

Our goal is to supply frontline workers, government employees as well as the general public with high-quality, American made masks at an affordable price.



125 Employees
Located in Houston

Purpose

Introduce BNX Converting, LLC (BNX) capabilities and solutions that will fulfill requirements, close capability gaps, and provide potential technological advancements in support of the COVID-19 Response Acquisition Task Force (DAF ACT) mission to provide relief, resilience, recovery, and stability to the nation in response to the COVID-19 pandemic.

Adam Ravat
CEO of BNX Converting LLC.

Rachel Leong
QC Lead
Mechanical Engineer

Guangli Ma
Lead Mechanical
Engineer

Yan Tong Ng Zhen
Quality Control Engineer

Alex Tan
Mechanical Engineer



Company Evolution

2010
Amcrest founded.
Importing and
distributing
security
cameras.



2014
Amcrest
established
manufacturing of
security cameras
in Asia in 2014.



2015
Acquired current 94,000
sqft manufacturing and
distribution facility in
Houston, TX

Aug, 2020
BNX Converting launched by
Amcrest in 2020 to establish US
based manufacturing of N95
masks.
Establishment of Duckbill
respirator line.



Sept, 2020
Acquired new 79,000 sqft
expansion manufacturing facility
in Houston, TX in (to meet
increased demand and surge
requirements)



Oct, 2020
Achieved KN95
certification on
Duckbill masks.



Oct, 2020
Completed NIOSH
Application and Issued
3-Digit Manufacturer
Number.



Jan, 2021
NIOSH N95 On-Site Audit
Completed Successfully and
Provisional Pass.

Sep, 2021
BNX becomes
Amazon's best seller
Made in USA N95
Respirator Mask.





**Medical
workers**



**Law
enforcement**



**Hospitality
workers**



**Teachers and
Educators**



**Construction
workers**

**To assist healthcare providers, first responders and
frontline workers in our nation's response to COVID-19**

#bnxmadeinusa

H95 **VERTICAL FOLDING RESPIRATOR** (HEADBAND)
95% PFE Filtration / 99% BFE Filtration




N95
 **MADE IN USA**
MADE IN USA FROM GLOBALLY SOURCED COMPONENTS

National Institute for Occupational Safety and Health
NIOSH
NIOSH-Approved - see insert



OEKO-TEX®
CONFIDENCE IN TEXTILES
STANDARD 100
21.HUS.74590 HOHENSTEIN HTTI


Vertical Folding Design


Blocks 95% Of Airborne Particles


Headband Attachment


Adjustable Metal Nosepiece

Multi-Layer Protection

- Spunbond Polypropylene
- Hot Air Cotton Polypropylene
- N95 Grade BFE99 Meltblown
- Hot Air Cotton Polypropylene
- Spunbond Polypropylene

#bnxmadeinusa

E95 **VERTICAL FOLDING RESPIRATOR** (EARLOOP)
95% PFE Filtration / 99% BFE Filtration




KN95
GB2626-2019

 **MADE IN USA**
MADE IN USA FROM GLOBALLY SOURCED COMPONENTS

OEKO-TEX®
CONFIDENCE IN TEXTILES
STANDARD 100
21.HUS.74590 HOHENSTEIN HTTI 


Vertical Folding Design


Blocks 95% Of Airborne Particles


Earloop Attachment


Adjustable Metal Nosepiece

Multi-Layer Protection

- Spunbond Polypropylene
- Hot Air Cotton Polypropylene
- N95 Grade BFE99 Meltblown
- Hot Air Cotton Polypropylene
- Spunbond Polypropylene

#bnxmadeinusa

F95 **FOLDING CUP STYLE RESPIRATOR**
95% PFE Filtration / 99% BFE Filtration




N95
 **MADE IN USA**
MADE IN USA FROM GLOBALLY SOURCED COMPONENTS

National Institute for Occupational Safety and Health
NIOSH
NIOSH-Approved - see insert



OEKO-TEX®
CONFIDENCE IN TEXTILES
STANDARD 100
21.HUS.74590 HOHENSTEIN HTTI


Folding Cup Style Design


Blocks 95% Of Airborne Particles


Headband Attachment


Adjustable Metal Nosepiece

Multi-Layer Protection

- Spunbond Polypropylene
- Hardshell Stiffer Polypropylene
- N95 Grade BFE99 Meltblown
- Spunbond Polypropylene

#bnxmadeinusa

C95 **CUP STYLE RESPIRATOR MASKS**
95% PFE Filtration / 99% BFE Filtration



N95
 **MADE IN USA**
MADE IN USA FROM GLOBALLY SOURCED COMPONENTS

National Institute for Occupational Safety and Health
NIOSH
NIOSH-Approved - see insert



Multi-Layer Protection

Spunbond Polypropylene
N95 Grade BFE99 Meltblown
Needlepunched Polypropylene



CUP Style Design



Blocks 95% Of Airborne Particles



Headband Attachment



Adjustable Metal Nosepiece

#bnxmadeinusa

A96-2 **HORIZONTAL FOLDING DUCKBILL RESPIRATOR**

95% PFE Filtration / 99% BFE Filtration




N95
 **MADE IN USA**
MADE IN USA FROM GLOBALLY SOURCED COMPONENTS

National Institute for Occupational Safety and Health
NIOSH
NIOSH-Approved - see insert



OEKO-TEX®
CONFIDENCE IN TEXTILES
STANDARD 100
21.HUS.74590 HOHENSTEIN HTTI


Horizontal Folding Design


Blocks 95% Of Airborne Particles


Headband Attachment


Adjustable Metal Nosepiece

Multi-Layer Protection
Spunbond Polypropylene
N95 Grade BFE99 Meltblown
Spunbond Polypropylene
Hydrophilic Spunbond Polypropylene



Quality, Safety, Innovation.

BNX's parent company, Amcrest, has been at the forefront of manufacturing security and safety products since its inception. BNX has leveraged this expertise and the necessary resources to fulfill its mission to manufacture the highest grade N95 NIOSH certified masks in order to support our great nation in response to the COVID-19 pandemic. BNX has a scalable and viable solution in place utilizing the highest grade materials, latest automated equipment as well as the talent, people, drive and unrelenting determination to fulfill our commitment.

01

**Competitively
priced** with
**highest level
quality**

02

Made in USA.
**Currently employ
125 people in
Houston, TX.**
Projected to add
at least 160 new
local jobs upon
expansion.

03

**Our focus is on
automated
manufacturing
in the US**

04

**Large investment in
automation of systems
and processes,** long
term sustainability and
viability in the US

05

**Wide array of products
(4 unique N95 mask
styles) to fit faces of all
shapes, sizes and
preferences.** We also have
Earloop KN95 vertical
folding and FFP2/CE
versions for convenience



06

Manufacturing capacity, we can scale up to 8million+ masks per month based on our current facility and machines

07

Supports local US jobs and building up strategic medical manufacturing facilities in Houston, TX

08

Profits are being used to re-invest in further US manufacturing. We have meltblown line on the way and are setting up for Air filter and water filter manufacturing in Houston.

09

Vertical integration for most competitive pricing and direct control over quality of input materials. We will be producing our own meltblown (the active filter ingredient) **in our masks and thus we will be able to pass those cost efficiencies on to our customers**

10

We have invested in 2 x TSI 8130A filter testing machines for our quality control. We have a third equivalent machine on the way

11

We are in the process of implementing our ISO13485 quality control system and will seek surgical 510K rating on all of our N95 masks

12

We have completed our NIOSH on-site Site Qualification visit and have been successfully granted provisional authorization for the production of our Duckbill N95 respirator

13

Operate out of a 94,000 sqft facility and acquired new 79,000 sqft adjacent facility to meet increased demand and surge requirements



Current Production Capacity & Scaling Requirements

October 2020
Establishment of
Duckbill Mask Line
with capacity up to
680,400
masks per month


December 2020
Establishment of
Cup Mask Line
with capacity up to
1,260,000
masks per month

January 2021
Establishment of
Vertical Headband Mask Line
with capacity up to
3,528,000
per month


January 2021
Establishment of
Folding Cup Mask Line
with capacity up to
1,296,000
masks per month



Projected Production Capacity



Testing Capacity





DEPARTMENT OF HEALTH & HUMAN SERVICES

NIOSH Reference: TN-24677
Mfr. Reference: BNX003

Centers for Disease Control
and Prevention (CDC)
National Institute for Occupational
Safety and Health (NIOSH)
National Personal Protective
Technology Laboratory (NPPTL)
626 Cochrans Mill Road
Pittsburgh, PA 15236-0070
Phone: 412-386-4000
Fax: 412-386-4051
April 23, 2021

Mr. Simon Ma
Chief Technology Officer
BNX Converting, LLC
16727 Park Row
Houston, TX 77084

Dear Mr. Ma:

The National Institute for Occupational Safety and Health (NIOSH) has reviewed the request referenced above and accepted March 12, 2021. This request is for approval of the H95W (white) and H95B (black) air-purifying filtering facepiece respirators for protections against particulates at a N95 filter efficiency level. The only difference between the respirators is the color. The complete respirator configurations are detailed on assembly matrix, file name *1113AMb.xls*, revision 2, dated: 04/08/2021.

This request is granted. Approvals are granted only for documentation written in the English language. It is the manufacturer's responsibility to correctly translate materials desired to be produced in languages other than English. Approval number **TC-84A-9315** has been assigned. This respirator is approved for protection against particulates at a N95 filter efficiency level.

The approval label is included as an attachment to this letter. The abbreviated labels have been accepted as submitted. The cautions and limitations are listed on the approval label. Only those assemblies affected by this request, or where new approval numbers are assigned, apply to this approval action. Production approval labels cannot include information on unapproved configurations.

The approved assemblies consist of the parts as listed on the approval label and the assembly matrix. Parts are to be marked with the numbers indicated on the approval label in a legible and permanent manner (marking cannot be removed without evidence of its previous presence).

This certificate of approval is not an endorsement of the respirator by NIOSH, and such endorsement shall not be stated or implied in advertisements or other publicity. However, you may publicize the fact that this respirator has met the requirements of Title 42, *Code of Federal Regulations*, Part 84 (42 CFR 84).

TC-84A-9315



DEPARTMENT OF HEALTH & HUMAN SERVICES

NIOSH Reference: TN-24896
Mfr. Reference: BNX004

Centers for Disease Control
and Prevention (CDC)
National Institute for Occupational
Safety and Health (NIOSH)
National Personal Protective
Technology Laboratory (NPPTL)
626 Cochrans Mill Road
Pittsburgh, PA 15236-0070
Phone: 412-386-4000
Fax: 412-386-4051
July 19, 2021

Simon Ma
Chief Technology Officer
BNX Converting, LLC
16727 Park Row
Houston, TX 77084

Dear Simon Ma:

The National Institute for Occupational Safety and Health (NIOSH) has reviewed your request accepted May 10, 2021. This request was for approval of the F95W and F95B N95 air-purifying filtering facepiece respirators. The respirators, identified by the color of the outer fabric, are part numbers F95W (white), and F95B (black). The only difference is the color. The complete respirator configurations are detailed on assembly matrix file name *1114AMa.xls*, revision A, dated: 4/21/2021.

This request is granted. Approvals are granted only for documentation written in the English language. It is the manufacturer's responsibility to correctly translate materials desired to be produced in languages other than English. Approval number **TC-84A-9362** has been assigned. This respirator is approved for protection against particulates at N95 filter efficiency level.

The approval label is included as an attachment to this letter. The abbreviated labels have been accepted as submitted. The cautions and limitations are listed on the approval label. Only those assemblies affected by this request, or where new approval numbers are assigned, apply to this approval action. Production approval labels cannot include information on unapproved configurations.

The approved assembly consists of the parts as listed on the approval label and the assembly matrix. Parts are to be marked with the numbers indicated on the approval label in a legible and permanent manner (marking cannot be removed without evidence of its previous presence).

This certificate of approval is not an endorsement of the respirator by NIOSH, and such endorsement shall not be stated or implied in advertisements or other publicity. However, you may publicize the fact that this respirator has met the requirements of Title 42, *Code of Federal Regulations*, Part 84 (42 CFR 84).

TC-84A-9362



FDA Registration #: 3017489417



Centers for Disease Control and Prevention (CDC)
National Institute for Occupational Safety and Health (NIOSH)
National Personal Protective Technology Laboratory (NPPTL)
626 Cochrans Mill Road
Pittsburgh, PA 15236-0070
Phone: 412-386-4000
Fax: 412-386-4051

The request is for approval of the model C95W (white) N95 air-purifying filtering facepiece respirator. The complete respirator configuration is detailed on assembly matrix file name *11154MB.xls*, revision B, dated: 06/21/2021.

The request is granted. Approvals are granted only for documentation written in the English language. It is the manufacturer's responsibility to correctly translate materials desired to be produced in languages other than English. Approval number TC-84A-9363 has been assigned. This respirator is approved for protection against particulates at N95 filter efficiency level.

The approved assembly consists of the parts as listed on the approval label and the assembly matrix. Parts are to be marked with the numbers indicated on the approval label in a legible and permanent manner (marking cannot be removed without evidence of its previous presence).

No changes may be made to any respirators and accompanying documentation without prior written approval of NIOSH. Requests for changes must be submitted to NIOSH and a modification of this approval must be granted before changes are made.

Enclosures



DEPARTMENT OF HEALTH & HUMAN SERVICES

NIOSH Reference: TN-24957
Mfr. Reference: BNX005

Simon Ma
Chief Technology Officer
BNX Converting, LLC
16727 Park Row
Houston, TX 77084

Dear Simon Ma:

The National Institute for Occupational Safety and Health (NIOSH) has reviewed your request accepted June 3, 2021. This request was for approval of the model C95W (white) N95 air-purifying filtering facepiece respirator. The complete respirator configuration is detailed on assembly matrix file name *11154MB.xls*, revision B, dated: 06/21/2021.

This request is granted. Approvals are granted only for documentation written in the English language. It is the manufacturer's responsibility to correctly translate materials desired to be produced in languages other than English. Approval number TC-84A-9363 has been assigned. This respirator is approved for protection against particulates at N95 filter efficiency level.

The approval label is included as an attachment to this letter. The abbreviated label has been accepted as submitted. The cautions and limitations are listed on the approval label. Only those assemblies affected by this request, or where new approval numbers are assigned, apply to this approval action. Production approval labels cannot include information on unapproved configurations.

The approved assembly consists of the parts as listed on the approval label and the assembly matrix. Parts are to be marked with the numbers indicated on the approval label in a legible and permanent manner (marking cannot be removed without evidence of its previous presence).

This certificate of approval is not an endorsement of the respirator by NIOSH, and such endorsement shall not be stated or implied in advertisements or other publicity. However, you may publicize the fact that this respirator has met the requirements of Title 42, *Code of Federal Regulations*, Part 84 (42 CFR 84).

No changes may be made to any respirators and accompanying documentation without prior written approval of NIOSH. Requests for changes must be submitted to NIOSH and a modification of this approval must be granted before changes are made.

Sincerely,

Jeffrey Peterson
Chief, Conformity Verification and
Standards Development Branch

Enclosures



DEPARTMENT OF HEALTH & HUMAN SERVICES

NIOSH Reference: TN-24584
Mfr. Reference: BNX002

Mr. Simon Ma
Chief Technology Officer
BNX Converting, LLC
16727 Park Row
Houston, TX 77084

Dear Mr. Ma:

The National Institute for Occupational Safety and Health (NIOSH) has reviewed the request referenced above and accepted February 4, 2021. This request is for approval of the A96-2 air-purifying filtering facepiece respirator for protections against particulates at a N95 filter efficiency level. The complete respirator configuration is detailed on the assembly matrix, file name *1112AMB.xls*, revision 2, dated: 04/07/2021.

This request is granted. Approvals are granted only for documentation written in the English language. It is the manufacturer's responsibility to correctly translate materials desired to be produced in languages other than English. Approval number TC-84A-9308 has been assigned. This respirator is approved for protection against particulates at a N95 filter efficiency level.

The approval label is included as an attachment to this letter. The abbreviated label has been accepted as submitted. The cautions and limitations are listed on the approval label. Only those assemblies affected by this request, or where new approval numbers are assigned, apply to this approval action. Production approval labels cannot include information on unapproved configurations.

The approved assembly consists of the parts as listed on the approval label and the assembly matrix. Parts are to be marked with the numbers indicated on the approval label in a legible and permanent manner (marking cannot be removed without evidence of its previous presence).

This certificate of approval is not an endorsement of the respirator by NIOSH, and such endorsement shall not be stated or implied in advertisements or other publicity. However, you may publicize the fact that this respirator has met the requirements of Title 42, *Code of Federal Regulations*, Part 84 (42 CFR 84).

No changes may be made to any respirators and accompanying documentation without prior written approval of NIOSH. Requests for changes must be submitted to NIOSH and a modification of this approval must be granted before changes are made.

Sincerely,

Jeffrey Peterson
Chief, Conformity Verification and
Standards Development Branch

Enclosures

TC-84A-9363

TC-84A-9308



FDA Registration #: 3017489417



Sponsor:
Adam Ravett
BNX Converting, LLC
16727 Park Row,
Houston, TX 77084

Sodium Chloride (NaCl) Aerosol Test Final Report

Test Article: A96
Study Number: 1332300-S01
Study Received Date: 18 Aug 2020
Testing Facility: Nelson Laboratories, LLC
6260 S. Redwood Rd.
Salt Lake City, UT 84123 U.S.A.
Test Procedure(s): Standard Test Protocol (STP) Number: STP0014 Rev 09
Deviation(s): None

Summary: This procedure was performed to evaluate particulate filter penetration as specified in 42 CFR Part 84 and TEB-APR-STP-0059 for requirements on a N95 respirator. Respirators were conditioned then tested for particle penetration against a polydispersed, sodium chloride (NaCl) particulate aerosol. The challenge aerosol was dried, neutralized, and passed through the test article at a concentration not exceeding 200 mg/m³. The initial airflow resistance and particle penetration for each respirator was determined.

According to 42 CFR Part 84.64, pretesting must be performed by all applicants as part of the application process with NIOSH. Results seen below are part of that pretesting and must be submitted to and accepted by NIOSH for respirator approval.

All test method acceptance criteria were met. Testing was performed in compliance with US FDA good manufacturing practice (GMP) regulations 21 CFR Parts 210, 211 and 820.



Trang Truong electronically approved for
Study Director
Curtis Gerow

24 Sep 2020 17:04 (+00:00)
Study Completion Date and Time



Study Number 1332300-S01
Sodium Chloride (NaCl) Aerosol Test Final Report

Results: The NIOSH N95 filter efficiency as stated in 42 CFR Part 84.181 is a minimum efficiency for each filter of ≥95% (≤5% penetration). The test articles submitted by the sponsor conform to the NIOSH N95 criteria for filter efficiency.

Test Article Number	Corrected ^a Airflow Resistance (mm H ₂ O)	Particle Penetration (%)	Filtration Efficiency (%)
1	17.1	0.779	99.221
2	18.1	1.52	98.48
3	17.3	3.16	96.84
4	17.2	0.629	99.371
5	17.3	1.76	98.24
6	17.4	2.04	97.96
7	17.2	3.61	96.39
8	17.3	1.73	98.27
9	17.6	1.69	98.31
10	17.9	0.418	99.582
11	17.7	0.584	99.416
12	19.1	0.490	99.510
13	18.3	0.528	99.472
14	17.5	0.799	99.201
15	17.5	0.383	99.617
16	18.4	0.581	99.419
17	17.8	0.380	99.620
18	18.1	0.298	99.702
19	18.6	0.478	99.522
20	17.5	0.898	98.102

^a The final airflow resistance value for each test article was determined by subtracting out the background resistance from the system.



Sponsor:
Adam Ravett
BNX Converting, LLC
16727 Park Row,
Houston, TX 77084

Determination of Inhalation and Exhalation Resistance for Air-Purifying Respirators Final Report

Test Article: A96
Study Number: 1331743-S01
Study Received Date: 15 Aug 2020
Testing Facility: Nelson Laboratories, LLC
6260 S. Redwood Rd.
Salt Lake City, UT 84123 U.S.A.
Test Procedure(s): Standard Test Protocol (STP) Number: STP0145 Rev 05
Deviation(s): None

Summary: This procedure was performed to evaluate the differential pressure of non-powered air-purifying particulate respirators in accordance with 42 CFR Part 84.180. The air exchange differential or breathability of respirators was measured for inhalation resistance using NIOSH procedure TEB-APR-STP-0007 and exhalation resistance with NIOSH procedure TEB-APR-STP-0003. The differential pressure technique is a simple application of a basic physical principle employing a water manometer differential upstream and downstream of the test material, at a constant flow rate.

According to 42 CFR Part 84.64, pretesting must be performed by all applicants as part of the application process with NIOSH. Results seen below are part of that pretesting and must be submitted to and accepted by NIOSH for respirator approval.

The inhalation resistance criteria as stated in 42 CFR Part 84.180 is an initial inhalation not exceeding 35 mm water column height pressure. The test articles submitted by the sponsor conform to this NIOSH criterion for airflow resistance.

The exhalation resistance criteria as stated in 42 CFR Part 84.180 is an initial exhalation not exceeding 25 mm water column height pressure. The test articles submitted by the sponsor conform to this NIOSH criterion for airflow resistance.

All test method acceptance criteria were met. Testing was performed in compliance with US FDA good manufacturing practice (GMP) regulations 21 CFR Parts 210, 211 and 820.

Results:

Test Article Number	Inhalation Resistance (mm H ₂ O)	Exhalation Resistance (mm H ₂ O)
1	12.4	13.7
2	13.8	13.0
3	13.7	13.6



Robert Dieker electronically approved for
Study Director
Curtis Gerow

04 Sep 2020 18:12 (+00:00)
Study Completion Date and Time

Filtration Report

Inhalation Exhalation Report



Certifications



FDA Registration #: 3017489417



Watch Video

#bnxmadeinusa

DEVICE REGISTRATION AND LISTING INFORMATION

Please review your registration.

Registration Date: Oct 16, 2020
Registration Status: Active, Waiting for Registration Number Assignment

Section 01		Type of Registration
Activities	MANUFACTURER COMPLAINT FILE ESTABLISHMENT	
Registration Number	Not Yet Assigned	
Section 02		Facility Name / Address Information
Name	BNX CONVERTING, LLC	
Address	BUILDING 1 16727 PARK ROW HOUSTON, TEXAS 77084 UNITED STATES	
Section 03		Owner/Operator Information
Owner/Operator Number	10078281	
Contact Name	ILENIA TAMAYO	
Company	BNX CONVERTING, LLC	
Address	16727 PARK ROW BUILDING 1 HOUSTON, TEXAS 77084 UNITED STATES	
Phone Number	(713) 936-2726	
Fax Number		
Email Address	Compliance@bnx.com	
Section 04		Official Correspondent Information
OC Contact Name	ILENIA TAMAYO	
OC Business Name	BNX CONVERTING, LLC	
Address	16727 PARK ROW BUILDING 1 HOUSTON, TEXAS 77084 UNITED STATES	
Phone Number	(713) 936-2726	
Fax Number		
Email Address	Compliance@bnx.com	
Section 05		Trade Names
Section 06		Listings
Listing #1		
Listing Number	D422096	
Listing Status	ACTIVE	
Created Date	Oct 16, 2020 12:18:03 PM ET	
Submission Number		
Proprietary Names	A96 Protective Mask	
Section 07		PIN/PCN Details
PIN	50304562	
PCN	21459164	

FDA Establishment Registration



Testing Report		Security website: www.fcl-sz.org.cn Security code: 6245716546	
Report No: ZFLJ2642376A		Page 1 of 8	
Applicant Information			
Applicant Name	: BNX CONVERTING,LLC.		
Applicant Address	: 16727 Park Row Houston TX 77084 USA		
Manufacturer	: BNX CONVERTING,LLC.		
Sample Information			
Sample Description	: A96 Duckbill mask		
Style No	: A96		
Brand	: BNX Converting LLC		
Sample Quantity	: 65 pieces		
Size specification	: 240mmX95mm		
Material Components	: 6 layers 1. 60gsm spunbond pp 2. 25gsm meltblown 3. 25gsm meltblown 4.15gsm spunbond pp 5. 25gsm skin friendly (hydrophilic) nonwoven 6. 25gsm skin friendly (hydrophilic) nonwoven		
- Sample Receiving Date : 2020-08-14			
- Report Date : 2020-08-26			
- The original sample is stuck on the last paper.			
Test Performed			
Judgement according to: GB 2626-2019 Respiratory protection—Non-powered air-purifying particle respirator			
- Selected test(s) as requested by applicants. For details, refer to attached page(s).			
Pronounce			
The results shown in this report refer only to sample(s) tested unless otherwise stated.			
All the items are performed in the standard conditions, except the noted cases.			
Except for the requirement of the client, the test results and the conformity judgement of this report do not take the uncertainty of the test results into account.			
Signed for and on behalf of CNTAC Testing Service Co.,Ltd.(Foshan)			
Approved by 			
广东省佛山市南海区西樵镇南方技术创新中心一、三、四楼1,364F, Nanfang Technology Innovation Center, Xiqiao Nanhai District, Foshan City, Guangdong, China Tel: (86 757)86850633/86806656 Fax: 86850633 Http://www.fcl.org.cn 若对检测报告有异议,应于收到报告之日起15日内向检测单位提出,逾期不予受理。			

CNTAC KN95 Test Report





© 2022 BNX Converting LLC. All Rights Reserved.
16727 Park Row Houston, TX 77084-5020
Phone: +1 713 936 2726 / +1 315 790 0466
Email: sales@bnx.com

For more information go to **www.BNX.com**

